



**QUOTATION**

DATE: August, 2008  
TO: SAMPLE QUOTATION

**Now with 5-Year Warranty**

REF: Vista VO2 Lab

We are pleased to offer the following Quotation #MASTER

ITEM	DESCRIPTION	UM	PRICE
<b>17605</b>	<b>Vista VO2 Lab Metabolic Measurement System</b>	<b>ea.</b>	<b>\$6,345.00</b>

Includes the following:

VISTA VO2 Lab with USB Interface, TurboFit-Windows Software for resting (REE) and Exercise Stress Testing, O2 Analyzer, automatic barometric pressure and temperature sensors; turbine flow sensor & mixing chamber. Also includes high-flow breathing valve for VO2max testing and air cushion mask starter kit. (Kit includes 20 masks in mixed sizes)

**Recommended Accessories**

1092	Calibration Syringe, 3-liter	325.00
16934 (*)	Calibration gas tank 16.6% O2 for approx. 20 calibrations	67.00
16939	Tank pressure regulator for above SD-size tank	68.00

**Total Recommended Accessories \$ 460.00**

**Recommended Accessories for Exercise Testing**

17038A	Polar Heart Rate Watch Interface (Imports heart rate into software, be sure to order P-1050 Polar watch, see below, if you do not already have a Polar watch)	295.00
P-1050	Polar FS1	59.95
VKIT-R	Reusable V-Mask Set: small, medium & large	365.00
R9066	Exercise Mouthpiece with Saliva Trap	37.00

**Total Accessories for Exercise Testing \$ 756.95**

**Other Options & Accessories**

16914 (*)	Calibration gas ED-size, 4% CO2, 16% O2, bal N2 (ED-size tank provides 5+ years supply)	159.00
16949	Tank pressure regulator for ED-size tank	175.00
17614	Mobile mini-cart for VO2 Lab, PC, printer & gas	750.00
17043	Digital Humidity Indicator	65.00
99999	On-site Installation/Instruction visit within 60 days of ship date at company convenience. (US only)	600.00
13620	Institutional PC-controlled treadmill, up to 400 pound load	6,345.00
OR		
17960	Corival PC-controlled ergometer to 999 watts	6,845.00
13610-6	Overhead suspension kit for breathing valve	195.00

**Replacement Resting Test Setups**

KM300	Starter Kit for disposable mask test setups, contains 3 samples of each below, headgear and 22mm tube adapter	149.00
KM301	Pack of 10 disposable test setups for resting metabolic measurements Each test setup includes a non-rebreathing valve, <u>LARGE</u> air cushion mask and connecting tube.	99.00
KM302	Pack of 10 disposable test setups for resting metabolic measurements Each test setup includes a non-rebreathing valve, <u>MEDIUM</u> air cushion mask and connecting tube.	99.00
KM302	Pack of 10 disposable test setups for resting metabolic measurements Each test setup includes a non-rebreathing valve, <u>SMALL</u> air cushion mask and connecting tube.	99.00

**Replacement Air-Cushion Masks**

KM201	Air cushion mask large adult (20/box)	\$89.95
KM202	Air cushion mask medium adult (20/box)	89.95
KM203	Air cushion mask small adult (20/box)	89.95
KM800	Head support for KM mask (spare)	25.00
KM900	Strap attachment plate (spare)	25.00

**Permanent (reusable) Face Masks**

Y2700	High-flow T-valve (spare)	260.00
R669170	"V-Mask" permanent silicone large face piece only	75.00
R669171	"V-Mask" permanent silicone medium face piece only	75.00
R669172	"V-Mask" permanent silicone small face piece only	75.00
R669173	"V-Mask" permanent silicone extra small face piece only	75.00
R200583	Adapter: V-Mask face piece to high-flow T-valve	95.00
R201418	Large Quick-Release headgear	69.00

Page 3 of 3  
Quotation #MASTER  
Aug 2008

Note: For high-volume testing we recommend you purchase additional valves, adapters & headgear.

PRICE INCLUDES FACTORY SERVICE

THE ABOVE PRICES ARE IN U S DOLLARS

Shipping schedule after receipt of order: 1-5 weeks depending on options purchased!

Terms: Net 30 days

Terms for Rental: See Rental Agreement

FOB: Ventura, CA

Valid: 60 days

Vacumetrics, Inc., VacuMed Division

By \_\_\_\_\_  
John Hoppe, President

\*17605.doc

P.S. Trade-in/Upgrade

Should you decide to upgrade to a full-featured VacuMed Metabolic Cart, we offer the following trade-in credit:

During the first 360 days from date of shipment: 80% of amount paid

During the next 360 days from date of shipment: 70% of amount paid