

Can your Metabolic Cart do This?

- ✓ **Metabolism Studies:**
Superior accuracy, results you can trust!
- ✓ **New Variables:** *VO₂sum, HRsum, Kcalsum, ΔHR, ΔW, ΔVE, True O₂, BPdia, BPsyst, BPmean, BPPP, RPP (double product), Watts calculation for all treadmill protocols. Plot against time or any other variable.*
- ✓ **Control Exercise Devices:**
Automatically adjust workload from PC
- ✓ **New RQ Check:** *Verifies proper gas analyzer calibration and software calculations during test independent of ventilation.*
- ✓ **Change Inspired Gas Concentration**
before or DURING test. (Think Altitude Training! A Low-O₂ delivery system is also available)
- ✓ **Customizable Screen Display:** *Create 6 different screen layouts and select at will during test.*
- ✓ **Tidal Breathing Display:** *Live on-screen and replay anytime after test.*
- ✓ **Expanded Data Averaging:** *All data may be averaged (filtered) breath-by-breath or at timed intervals. Example: Display a 30-second average every 10 seconds.*
- ✓ **Hill Protocol:** *Customize height and distance between hills.*
- ✓ **Linear Treadmill Ramp Protocol** *converts ramp speed and elevation to watts.*
- ✓ **Test Overlay:** *Pick tests from same or different subjects and overlay them graphically. (Think Before and After)*
- ✓ **New Fitness Analysis Report** *includes bar graphs, personalized lactate threshold-based training heart rate range and 4-week Training Schedule.*
- ✓ **Print Explanation** *for Exercise or Resting test so that test subject understands the results.*
- ✓ **Customizable Reports:** *Create groups of pages that will print with clicking a single button.*
- ✓ *Many more features to impress you and your customers, ask for demo disk.*